

### **BS-2044 Biological Microscope**



BS-2044B BS-2044T

#### Introduction

BS-2044 series microscopes are high quality biological microscopes, which are specially designed for biological and medical research and teaching experiments for colleges, universities, laboratories and related institutes. With infinity color correction optical system and excellent Koehler illumination system, BS-2044 can get uniform illumination, clear and bright images at any magnification. These microscopes can be used for teaching experiments, pathological examinations and clinical diagnosis. With outstanding functions, excellent cost performance, easy and comfortable operation, BS-2044 series microscopes present over expected and splendid micro images.

#### **Feature**

- 1. Infinite Color Corrected Optical System provides sharp and comfortable images.
- 2. Wide-field high eye-point eyepieces and plan achromatic objectives make the effect of fluorescence observation more perfect.
- 3. Outstanding appearance design and ergonomics structure design, long time use without feeling tired.
- 4. With the design of security lock and safe limits, it is much more secure, stable and can be kept for long.
- 5. Various microscopic examinations can be fulfilled, such as brightfield, dark-field, phase contrast, fluorescence, simple polarizing and so on.
- 6. LED fluorescence excitation illumination, breaking through the traditional type, is more stable, low radiation and long working life. Special fluorescent filters for tuberculosis examination is available.

#### **Application**

BS-2044 series biological microscopes are ideal instruments for biological, pathological, histological, bacterial, immune, pharmacological and genetic fields. They can be widely used in education, medical and sanitary establishments, such as schools, hospitals, clinics, laboratories, medical academies, colleges, universities and related teaching labs and research centers.



# Specification

Item		Specification	BS-2044B	BS-2044T
Optical System	Infinite Color Corrected Optical	System, parfocal distance 45mm	•	•
Viewing Head	Seidentopf binocular head, 30°	inclined, Interpupillary 50-75mm, 360° rotatable,	_	_
	eyepiece tube: Φ30mm		•	0
	Seidentopf trinocular head, 30° inclined, Interpupillary 50-75mm, 360° rotatable,		0	•
	Fixed light splitting ratio: Eyepie	Fixed light splitting ratio: Eyepiece:Trinocular=8:2, eyepiece tube: Φ30mm		
	Seidentopf trinocular head(ded	Seidentopf trinocular head(dedicated for fluorescence), 30° inclined,		
	Interpupillary 50-75mm, 360° rotatable, Fixed light splitting ratio:		0	0
	Eyepiece:Trinocular=5:5, eyepiece tube: Φ30mm			
Eyepiece	High eyepoint wide field plan eyepiece PL 10×/22mm with adjustable diopter ±5		•	•
	High eyepoint wide field plan eyepiece PL 10×/22mm with adjustable diopter ±5,		0	0
	with eyepiece micrometer			
	Eyepiece pointer		0	0
	Eyepiece micrometer		0	0
		2×, NA=0.06, WD=5.0mm (should work with	0	0
		frosted filter)		
		4×, NA=0.10, WD=15.0mm	•	•
	Infinite Plan Achromatic	10×, NA=0.25, WD=10.8mm	•	•
	Objectives	20×, NA=0.45, WD=1.5mm	0	0
		40×(S), NA=0.65, WD=0.8mm	•	•
		60×(S), NA=0.85, WD=0.3mm	0	0
		100×(S, Oil), NA=1.25, WD=0.21mm	•	•
Objective		10×, NA=0.25, WD=12.1mm	0	0
	Infinite Plan Phase Contrast	20×, NA=0.40, WD=1.56mm	0	0
	Objective	40×(S), NA=0.65, WD=0.36mm	0	0
		100×(S, Oil), NA=1.25, WD=0.21mm	0	0
		4×, NA=0.13, WD=16.4mm	0	0
	Infinite Plan	10×, NA=0.30, WD=8.1mm	0	0
	Semi-Apochromatic	20×, NA=0.50, WD=2.0mm	0	0
	Fluorescence Objectives	40×(S), NA=0.75, WD=0.74mm	0	0
		100×(S, Oil), NA=1.28, WD=0.14mm	0	0
Nosepiece	Reversed Quadruple Nosepiece		•	•
	Reversed Quintuple Nosepiece		0	0
Stage	Double Layers Mechanical Stage 150mm×140mm, moving range 76mm×50mm,			
	double slide holder, accuracy: 0.1mm		•	•
	Rackless Double Layers Mechanical Stage 150mm×162mm, moving range			
	76mm×50mm, double slide holder, accuracy: 0.1mm, Wear-resistant and		0	0
	anti-corrosion treatment on the stage surface			
Condenser	N.A.1.25 Koehler illumination condenser (with plug-in phase contrast and dark		•	•
	field plate slot), condenser preset center and height adjustable			
	N.A.1.25 Kohler Illumination Turntable Phase Contrast Condenser Set for 10X,		0	0
	20X, 40X, 100X Phase Contrast Observation			
	N.A.1.25 Kohler Illumination Tu	rntable Phase Contrast + Dark Field Condenser	0	0




# **BestScope International Limited**

	bestscope inter		
	Set for 10X, 20X, 40X, 100X Phase Contrast and 4X-40X Dark Field Observation		
Focusing	Low position coaxial focusing system, moving range 30mm, with upper limit and tightness adjustment, fine division 0.002mm		
			•
	Adaptive 100V-240V, AC50/60Hz wide range voltage, single high brightness 3W		
Transmitted Illumination	LED lamp (preset center, Cold/warm color temperature optional, cold color		
	temperature 4750K-5500K, warm color temperature 2850K-3250K), light intensity can be continuously adjusted		•
Mercury			
Reflected	mercury bulb (OSRAM/ Chinese brand)		0
Illumination			
	B1 band-pass type fluorescence module, with intensity adjusting knob, and	0	
	switch knob for bright field and fluorescence, central wavelength:470mm		0
	G1 band-pass type LED fluorescence module, with intensity adjusting knob, and		
	switch knob for bright field and fluorescence, central wavelength: 560mm	0	0
	B2 long-pass type fluorescence module, with intensity adjusting knob, and		
	switch knob for bright field and fluorescence, central wavelength:470mm	0	0
LED Fluorescent	UV2 ultraviolet long-pass type LED module, with intensity adjusting knob, and		0
Reflected	switch knob for bright field and fluorescence, central wavelength: 385mm	0	
Illumination	UV4 ultraviolet long-pass type LED module, with intensity adjusting knob, and		0
	switch knob for bright field and fluorescence, central wavelength: 365mm	0	
	B4 LED fluorescence module dedicated for TB, with intensity adjusting knob, and switch knob for bright field and fluorescence, central wavelength: 455mm		0
	Other various LED module for option, which can be custom-tailor according to the needs of clinical diagnosis.		
			0
Filters	Blue filter Φ45mm	0	0
	Green filter Φ45mm	0	0
	Yellow filter Φ45mm	0	0
	Neutral filter Φ45mm	0	0
Polarizing Kit	Polarizer	0	0
	Analyzer	0	0
Dark Field Plate	Dark field insert plate (used for 4×-40× objectives)	0	0
Centering		_	_
Telescope	Centering TelescopeΦ23.2mm (used with phase contrast plate and objective)	0	0
Phase Contact Plate	10×, 40× Phase Contact Insert Plate (used for 10×, 40× phase contrast objectives)	0	0
	20×, 100× Phase Contact Insert Plate (used for 20×, 100× phase contrast		0
	objectives)		
C-mount Adapter	0.35× C-mount adapter, adjustable	0	0
	0.5× C-mount adapter, adjustable		0
	1× C-mount adapter, adjustable	0	0
	Trinocular tube for digital eyepiece (Φ23.2mm)	0	0
Packing	1 set/carton, 58x56x28cm, GW: 10kgs, NW: 8kgs	•	•


Note: ● Standard Outfit, ○ Optional



## Sample Images



### **Dimension**



Unit:mm